

MIXPAC™

Two-Component System

18 to 75mL



FILLING GUIDELINE
Version 1.1



DENTAL

Contents

		Page
1	General	3
	1.1 Scope of Content	3
	1.2 System Name	3
	1.3 System Provider	3
	1.4 Terminology	4
	1.5 Symbols	4
2	System Overview	5
3	Basic Cartridge Filling Techniques	6
4	Quality of the Fill	7
5	Piston Insertion	8
6	Common Filling Problems	9
7	Front Filling	11
8	Filling Volume	12
9	Large Air Inclusions	13

1 General

1.1 SCOPE OF CONTENT

This document is intended as a general guideline for the process of filling and piston insertion. The information within provides a solid foundation to achieve high quality results. The guideline is applicable to manual operation as well as automated filling lines but the result of the filling process must in either case be verified by the customer (storage and application).

medmix provides the system. The receiving party shall be responsible and liable for the usage of the system and the materials applied therewith .

NOTE



In case of system compatibility issues, please do not hesitate to contact your local medmix Sales Manager.

We are grateful for any improvements suggestions and all constructive comments.

1.2 SYSTEM NAME

MIXPAC™ Two-Component System 18 to 75mL

1.3 SYSTEM PROVIDER

medmix Switzerland AG

Ruetistrasse 7

9469 Haag

Switzerland

Tel.: +41 81 414 70 00

E-Mail: mixpac@medmix.com

www.medmix.swiss/dental

1.4 TERMINOLOGY

The following terms are used in this filling guideline:

Multicomponent Cartridge Systems; Cartridge, Plunger, Pistons, Cap

Piston setting equipment; Air venting profile (AVP) aka venting shim

Dispensing (discharge material from the cartridge through the mixing tip– using a dispenser)

1.5 SYMBOLS

Please observe the meaning of the following symbols:

CAUTION



Topic

Please be aware of possible system failure if instructions are not followed.

NOTE



Additional information for specific application

2 System Overview

The application system consists of a multicomponent cartridge with preinstalled cap and matching pistons.



1. Cap
2. Separated outlets
3. Cartridge
4. Pistons with o-ring
5. Pistons with lip and o-ring

NOTE



Please refer to the MIXPAC™ Dental Products catalogue for detailed information.

3 Basic Cartridge Filling Techniques


Air trapped within the cartridge system can drastically degrade the performance of the two component material. The air inhibits the ability of the cartridge to meter the A- and B-side materials accurately because the air behaves like a shock absorber leading to a few problems; lead lag issues, poor mixing in the static mixer, and cross contamination.

The difficulty inherent with filling a cartridge air free can be mitigated by using sound techniques and properly designed equipment. Some basic techniques and/or ideas are shown below.

CAUTION



Air entrapment
Air trapped within the cartridge system can affect the performance of the two component material.



1. A properly made filling mandrel fits tightly to the inner diameter of the cartridge and provides the most effective means of filling. It allows a small amount of air to escape from the side. The material is pumped through the mandrel and the cartridge is pushed off the mandrel. In this scenario it is easy to add a limit switch to control the filling height. Additionally, it takes the control away from the operator as you will see in scenarios below.
2. Using a standard size filling nozzle is OK but is not as reliable as the situation found using a filling mandrel. The speed at which the filling nozzle is removed must be timed with the flow rate to minimize the risk of air entrapment. This setup works best if one is using a premeasured shot based system to fill the cartridge.
3. In this scenario the filling nozzle is moving faster than the material is deposited. It often leads to a large air bubble trapped towards the wall of the cartridge tube. With lower viscosity materials this scenario will create many very small bubbles as the fluid entering the chamber is very turbulent and entraps air.
4. In this scenario the filling nozzle is moving slower than the material is deposited. The material engulfs the end of the filling nozzle and creates a void as the nozzle is pulled from the cartridge. It often leads to large air bubbles being trapped at the center of the tube.

4 Quality of the Fill


It is very important that the A and B materials are filled to equal levels to insure the proper ratio can be achieved. If the materials are filled to different heights it drastically increases the likelihood of an off-ratio condition, cross contamination, disappointment from an end user, and possibly environmental waste. Below are some examples of some filled cartridges ranging from ideal to poor. These examples are intended to be used as a general guideline to achieving high results with your filling challenges.

CAUTION



Filling Level Equilibrium

Uneven filling levels and or filling surface might lead to air entrapment at the piston.



1. Both of the filling levels in the first scenario are equal, there is no air trapped within the material and the profile at the top of both A- and B-sides filling is slightly domed upwards.
2. In this scenario it can be seen that the filling heights are not equal as there is more blue material than yellow. This indicates that the process is neither accurate, nor repeatable. It will create a risk for cross contamination and lead lag problems at the start of dispensing the cartridge. However, there are no air-bubbles trapped in the material and is therefore rated OK.
3. In this scenario the filling is littered with small air bubbles that may not collect (dependent on filling viscosity). The presence of the air bubbles will cause metering problems. The filling profile for the yellow side is very poor and material at the upper left edge of the cartridge could inhibit the sealing function of the piston lips. Users are advised to avoid this.
4. In this scenario there is a large air bubble trapped on the B-side which will cause metering problems. Also, the profile of the B-side filling is concave and might inhibit the pistons ability to vent all of the air. Again, this type of scenario will produce poor metering results and should be avoided.

NOTE



See also " Effects of air inclusions - Multicomponent cartridge system"

5 Piston Insertion

18 to 75 mL cartridge systems can be filled from the front as well from the back. Influencing factors are piston design, filling media viscosity as well as process capabilities/ degree of automation. In each case it must be ensured to fully evacuate the air trapped between the piston and the filling.


The sections below lists the factors to be considered for an air-free filling process, without contaminating and/or damaging the sealing features. The illustrations show a generic multicomponent system 50 mL cartridge system.

CAUTION



SEALING FEATURES

Damage to or contamination of the sealing features will subsequently lead to leakage during transportation or storage of the filled cartridges and/or system failure upon usage.




1. Install cartridge into a properly designed guide block and then place piston into said block so that its rib structure can be seen from the top. The piston should fit down into the block and be square with the cartridge.
2. The guide block will press down onto the cartridge clamping it into place. Additionally, a venting shim will be guided into the cartridge.
3. Push the insertion plunger down forcing the piston into the cartridge body.
4. Carefully keep pushing the plunger until all of the air is evacuated. No material should leak along the shim.
5. The Insertion Plunger returns back to the home position. The guide block also returns and pulls the venting shim from between the cartridge and the piston.

Some highly automated equipment use vacuum technology to evacuate the air below the piston. This makes the use of a shim redundant and eases the process. Please refer to "Manufacturers of filling equipment" to identify an experienced partner in your area.


6 Common Piston Insertion Problems

Below are some further examples of common filling scenarios that will be encountered in any organization. Even further below there are explanations of each scenario and how to possibly rectify and/or avoid the situation as they are not ideal.




Leaking Cartridge

Never try and insert the piston at an angle as it can lead to damage to the sealing features. If the O-ring becomes damaged the cartridge will likely leak either in storage or during dispensing. It is always best to use a guide block and/or setting device.



Piston not functioning correctly


The piston in this scenario is upside down. The piston will not function correctly if oriented in this direction. The correct orientation can be identified by checking for the pistons rib structure.



O-ring disengages from original groove


The O-ring has rolled out of the original groove during insertion. Normally this will result in a leak. The piston should be removed and replaced with a new one.

Air is not able to escape




This scenario shows that the piston has been inserted into the cartridge such that not all air was allowed to escape from underneath the piston. This cartridge will likely not meter accurately leading to poor performance of the two component material. See also "Multicomponent systems - Effects of air inclusion"

Overfilling cartridge



Do not overfill the cartridge as it may leak and/or not fit into the dispensing tool.

Material leakage beyond the O-ring seal




In this scenario the plunger pushes too long and/or with too much force on the material such that it pushes up beyond the O-ring seal. Some fillings will tolerate small amounts of material getting past the seal but others will not.

7 Front Filling

medmix does not promote nose filling but it can be an applicable process under certain circumstances for cartridges sizes up to 50 mL. Below is a representation of a well-designed nose filling process.

CAUTION


Piston tilt
Keep the piston square to the cartridge when filling to prevent leakage



- 1) Filling Nozzle
- 2) Cartridge
- 3) Guide Plate
- 4) Pistons
- 5) Plunger
- 6) Air Cylinder

1. The cartridge is held securely in place with an appropriate alignment fixture. The outlets are not yet connected to the filling nozzles unless a venting/vacuum feature is available in the nozzle. The plungers are aligned with the cartridge and the setting process is initiated. Please ensure that the pistons are arranged concentrically with the openings. Please also ensure to select the correct type of piston and that the sealing features are not pre-damaged
2. The pistons are pushed towards the front of the cartridge tubes. The necessary force is dependent on the tube diameter, as well as the cartridge and piston material combination. The placing force and speed must be well controlled and verified by the customer in order to prevent damage to the sealing features (note: F1 up to 80N). After the air has been purged from the tubes, the cartridge outlets are connected to the filling nozzle. Some air will be entrapped at the interface between nozzle and pistons. For air-free filling, a venting or even better a vacuum mechanism should be integrated into the filling nozzle.
3. The filling nozzle is activated. The materials are forced into the A- and B-side tubes of the cartridge. The material pushes the plungers towards the back of the cartridge. To optimize the filling process, applying a light counterforce on the plunger may be beneficial. (note: F2 ~20N).
4. The filling volume must not exceed the net volume specified. The material flow has to be stopped with a dosing pump or limit switch. Using a mechanical end stop for the plunger is not recommended in order to prevent high compressive load in the system. An integrated valve function in the nozzle prevents material overflow and contamination of the sealing features.

8 Filling Volume



- AL** Assembly length
- CL** Cartridge length
- CW** Cartridge width
- ID-A** Inner diameter A-side
- ID-B** Inner diameter B-side
- PD** Piston depth

Size [ml]	Ratio	CL [mm]	AL [mm]	CW [mm]	ID-A [mm]	ID-B [mm]	PD [mm]
18	1:1	46	66.2	44.6	19.1	19.1	≥3
25	1:1	78	99	37.5	15.8	15.8	
30	1:1	67	87.6	44.6	19.1	19.1	
50	1:1	102	122.6	44.8	19.1	19.1	
	2:1	102	104	44.9	22.2	15.7	
	4:1	102	112	43.6	24.2	12.1	
	10:1	103	122.6	43.6	13.3	6.65	
75	1:1	148	168	44.8	19.1	19.1	

source: internal data

The dimensional information contained within this table is meant as a guideline. Valid information for the specific system in use (system drawing) can be requested from your local medmix Sales Manager.

9 Large Air Inclusions

Air trapped within the cartridge system can drastically reduce the performance of a two component material. The air inhibits the ability of the cartridge to deliver an accurate ratio between the A- and B-side materials. As the air behaves like a shock absorber, it can cause one or several of the following problems; lead lag issues, poor mixing, cross contamination and hence, potential degradation of the two component material.

Please refer to the separate info sheet for an illustrated description of the mechanics. The customer has to verify the amount of air entrapped in the system by his filling process as well as the effects of the latter on the mixing quality.

10 Technical Data

The technical data needed for your cartridge system might differ due to filling material specific needs. Please contact your local medmix sales manager in order to provide you the necessary information for your application.

11 Available Tools and Drawings from medmix

In order to vent pistons medmix offers venting shims (AVP's) for various cartridge diameters. Also basic piston setting tools (manual and pneumatic) are available from partners. Please contact your local medmix MIXPAC sales manager in order to provide you the necessary information for your application.

12 Applicable Documents

Effects of air inclusions - Multicomponent cartridge system
Manufacturers of Filling Equipment - Dental Products

EUROPE

medmix Switzerland AG

Ruetistrasse 7
9469 Haag
Switzerland

Tel.: +41 81 414 70 00
mixpac@medmix.com
www.medmix.swiss

USA, SOUTH AMERICA, CANADA AND MEXICO

medmix US Inc.

8 Willow St.
Salem, NH 03079
USA

phone +1 603 893 2727
mixpac@medmix.com
www.medmix.swiss

ASIA PACIFIC

medmix China Ltd.

68 Bei Dou Lu, Minhang
Shanghai 200245
China

phone +86 21 6430 68 68
mixpac@medmix.com
www.medmix.swiss



Swiss Quality
at the Highest Level



Legal Notice: Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC, T-Mixer, and Colibri are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland 2022