MIXPAC™






DSA / DS Dispensers

DENTAL

DSA Dispenser

DS Dispenser

en

1. Intended purpose
This non-sterile, manual, hand-held device is intended for dispensing two-component dental materials (e.g. elastomeric impression materials, dental cements) intra- and extraorally from a MIXPAC™ Cartridge and MIXPAC™ Mixing tip. This is a reusable device.

2. Indications for use
Indicated for use with dental materials such as elastomeric impression materials and dental cements (e.g. resin-based cement, Bis-Acrylic resin, resin-modified glass ionomer cement).

3. Contraindications
There are no known contraindications.

4. Intended Patient Population
The definitive patient group is diagnosed on a case-by-case basis by a licensed dental professional.

5. General Information
This dispenser is delivered non-sterile.

6. Using the product

6.1 Preparation

6.2 Application

7. Reprocessing

7.1 Disassembly

7.2 Limitations on processing
Any dispenser showing visible damage or wear must be exchanged.

End of life of a reusable dispenser is normally determined by wear and damage due to use.

Use-life of the dispenser is exceeded if visual inspection shows unacceptable deterioration like wear, discoloration, cracks, etc.

7.3 Initial treatment at the point of use

To prevent any contamination from drying it is recommended that reprocessing takes place not later than 2 hours after use.

Keep the dispenser moist after use to prevent soil from drying on them.

Contaminated dispenser should be transported to the area for cleaning in a clean and closed box/container that avoids spreading of the contamination.

7.4 Manual cleaning

Use a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) to clean all surfaces of the plunger for 30 seconds. Replace wipe, as necessary, if it becomes soiled.

In order to remove residue on the plunger, use a soft bristle brush or a toothbrush and brush for 30 seconds. If necessary, increase time until no visible soil is observed.

Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces of the plunger two (2) times to remove residual detergent agent. Use additional wipes as needed.

Place a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) through the gaps of the plunger and move the wipe back and forth for 30 seconds in the gaps to clean the plunger. Replace wipe, as necessary, if it becomes soiled.

In order to remove residue in the gaps of the plunger, use a soft bristle brush or a toothbrush and brush for 30 seconds. If necessary, increase time until no visible soil is observed.

Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces of the plunger two (2) times to remove residual detergent agent. Use additional wipes as needed.

Place a clean, lint-free wipe, dampened with a cleaning agent containing isopropyl alcohol (e.g. CaviWipes 17%) through the gaps of the dispenser and move the wipe back and forth for 30 seconds to clean the gaps. Replace wipe, as necessary, if it becomes soiled.

In order to remove residue in the grooves of the dispenser, use a soft bristle brush or a toothbrush and brush for 140 seconds. If necessary, increase time until no visible soil is observed.

Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces of the dispenser two (2) times to remove residual detergent agent. Use additional wipes as needed.

7.5 Manual disinfection

Use a clean, lint-free wipe, dampened with a disinfectant containing isopropyl alcohol (e.g. CaviWipes 17%) to allow all surfaces and gaps of the plunger and dispenser to remain visibly wet for a minimum of five (5) minutes. Use additional wipes to ensure all surfaces and gaps remain wet for the five (5) minute contact time.

Use a clean, lint-free wipe, dampened with purified water to wipe all surfaces and gaps of the plunger and dispenser for at least two (2) times to remove residual disinfectant. Use additional wipes as needed.

7.6 Drying

Allow the plunger and dispenser to dry until visibly dry.

7.7 Maintenance, Inspection and Testing

Inspect the dispenser and plunger under normal lighting:

- Any dispenser or plunger showing visible damage or wear must be exchanged.
- Repeat manual cleaning and inspection if the dispenser or plunger are not visibly clean.

7.8 Packaging

- Reassemble the dispenser and plunger.
- Place the dispenser in double sterilization packaging (pouches made of paper and plastic film), conforming to ISO 11607 and ANSI/AAMI ST79.
- For USA, use FDA-cleared pouches.
- Physical-chemical indicators may be used on or in the packaging system

7.9 Sterilization

Use the sterilization parameters provided in table below which have been validated:

Procedure	Dynamic-air-removal Steam Sterilization Cycle	
	USA	Europe
Exposure Time	4 minutes	Min. 3 minutes
Temperature	132 °C	134 °C
Drying Time	20 minutes	20 minutes

7.10 Storage

- Keep the dispenser in sterilization packaging, in a dry and clean environment
- Sterility cannot be guaranteed if packaging becomes open, damaged, or wet.
- Check the packaging and the dispenser before using it (packaging integrity, no excessive humidity and validity period).

8. Accessories

For the DS-Dispenser a Spacer Insert is available as an accessory.

8.1 Spacer Insert for DS-Dispenser

The spacer insert is used to reduce the distance of the handle and lever for better usability for users with small hands.

9. Product Disposal

Dispose according your local disposal requirements for dental devices.

If any serious incident has occurred, please report it to the legal manufacturer and the competent authority of the Member State in which the incident occurred.

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