



Imagine the Possibilities

FleXa™ is the ideal platform for product differentiation. medmix offers multiple tiers of **Customization & Branding**:



Off-the-shelf solution:

Choose from medmix's standard palette with over **150+ medical device-grade colors**.



Custom colored accessories:

Match your own colors for instant recognition, whether it's your corporate color, brand colors or identification colors for instance for tooth colors.



Bespoke syringe design:

Co-developed solutions tailored to your needs, including differentiated shapes, foil printing for a **unique premium look and feel**.

Quality, Innovation and Differentiation is in our DNA.

Accelerate Innovation with our end-to-end expertise.

At medmix, we partner closely with our customers to co-create solutions that elevate user experience. We accelerate innovation with comprehensive, end-to-end expertise.

Let's create your next market winner – Together!

Are you ready to make your next big idea move faster with medmix? Contact your sales representative now.

medmix Switzerland AG
Ruetistrasse 7
9469 Haag, Switzerland
P +41 81 414 70 00
www.medmix.swiss



MIXPAC and FleXa are trademarks of medmix Switzerland AG. Candy Colors is a trademark belonging to medmix Switzerland AG, and the Candy Colors do not mean diameter or size. The Candy Colors stand for medmix Switzerland AG as the source of the products, the MIXPAC brand, and Swiss quality. medmix is trademark of medmix AG. © medmix Switzerland AG 2026


The features and claimed performance presented are part of ongoing development and are subject to final validation. This information is preliminary, may change without prior notice, and does not constitute a binding commitment or guarantee of future availability or performance.



FleXa™

The syringe system that adapts to your brand and elevates every dental experience.

Make It Personal



000000-00_FleXa 01.26 en

What is FleXa™ ?

FleXa is a highly customizable, syringe system with advanced usability features, designed for a wide spectrum of one-component dental materials.

Universal syringe platform

FleXa is designed for premium dental materials. It is a silicone oil-free syringe system, developed for one-component gels, pastes, and liquids, where precision and control matter most for example:*



Fissure sealants, etchant gels and liquids



Flowable composites, veneer cements (LC), cavity liners and pulp capping materials



Endodontic dressings and internal bleaching pastes



In-office bleaching gels and gingival barrier/block out gel




Retraction paste and caulking paste



* Always needs OEM validation of product compatibility.

FleXa™ Syringes for single-component dental materials



Shared dimensions between volumes

allows for interchangeability between accessories

Dream-Team for composites and veneer cements:

Combine our Bubble-free dispensing tip with FleXa for enhanced precision.



Silicone oil free syringe

smooth application without silicone oil residues

Available in 4 sizes



0.75 mL 1 mL 2 mL 3 mL