



Twistpac™ Unit Dose

Superior Filling. Efficient Container.



Improving the dental workflow.

Clean. Reliable. Efficient

Designed to **enhance your product portfolio**, simplify workflows, and **deliver real value** to your customers.

Hygienic by design

Unit dose system minimizes the risk of cross-contamination.

User-friendly design

Easy to open due to predetermined breaking point: break, dip and apply.

Integrated pedestal for convenient working during dental treatment.

Enhanced efficiency

Accurate dosing for time-saving application.



For best results, pair your Twistpac™ Unit Dose with ZerofloX™ Micro Applicator

Your material, perfectly packaged

Experience the perfect synergy of innovative Unit Dose packaging and expert filling service.

Designed for liquids like **primer, adhesive, desensitizer, silver diamine fluoride, and similar solutions.**



Premium branded

with customized pad printing.

Precisely filled

up to 0.25 mL in every unit.

Securely sealed

by ultrasonic welding.

No additional foil wrapping required

- High quality Material COP (Cyclic Olefin Polymer)
- Excellent material protection from UV light
- Developed to protect product stability and integrity

Swiss-made Twistpac™ containers

for uncompromising accuracy.



Swiss-made quality.

High-speed filling lines in Germany

built for reliability and consistency.




German-engineered performance.

Twistpac™ Unit Dose From Vision to Market: Your Filling Service Roadmap.

1 Product Concept

Prior to the first trial run, we define all key product parameters. Together with your team, we coordinate **filling volume, unit dose color, printing layout and secondary packaging**—tailored to the requirements of your bulk material and market needs.



2 Trial Run

In the first project phase, we test your bulk material on our filling equipment to **determine optimal process parameters, define tolerances, and produce prototypes** for stability and market testing.

3 Serial Quality

Once trial samples are approved and specifications finalized, we begin **serial filling with a full qualification and validation process** to ensure consistent batch quality.



4 Finished Goods

As an added service, we help optimize your secondary packaging, source your selected chip-boxes and IFUs based on your artwork, and **assemble a complete 'ready to sell' kit**.

Expertise that elevates your product:

With over **20 years of experience** in filling a wide range of dental materials, our specialists help define the optimal process for your formulation — ensuring safety, stability, and perfect usability.

Tailored to your needs.



Use customized labeling

Individualized pad printing on every container.



Shine with your brand colours

Exclusive colour options possible starting from a minimum order quantity of 200,000 units.



Select your bagging solution

Bagged in a size optimized, resealable bag labeled with all the individual information you need (e.g. Logo, UDI, LOT, expiry Date).



Define your packaging needs

Individual box layout, including IFUs, technique cards, micro applicators.

Make it Yours.



medmix Switzerland AG | Ruetistrasse 7 | 9469 Haag, Switzerland | www.medmix.swiss
medmix Deutschland GmbH | Zum Kesselort 53 | 24149 Kiel, Germany

MIXPAC, Twistpac and ZerofloX are trademarks of medmix Switzerland AG. Candy Colors is a trademark belonging to medmix Switzerland AG, and the Candy Colors do not mean diameter or size. The Candy Colors stand for medmix Switzerland AG as the source of the products, the MIXPAC brand, and Swiss quality. medmix is trademark of medmix AG. © medmix Switzerland AG 2026

Learn more here.

