


Flexible Steel™ IRRIGATION NEEDLES



FLEXIBLE STEEL IRRIGATION NEEDLES

transcendent™ Flexible Steel Irrigation Needles are double-sided vent irrigation needles with Luer Lock that increase safety and efficacy of endodontic irrigation. The full working length can be easier achieved and effectively irrigated even in curved root canals.



ORDER NO.	DESCRIPTION	SIZE	COLOR	QTY/BOX
162066	23 G	0.6 x 25 mm	blue	40
162067	25 G	0.5 x 25 mm	orange	40
162068	27 G	0.4 x 25 mm	grey	40
162069	30 G	0.3 x 25 mm	yellow	40

PRODUCT FEATURES AND BENEFITS

Closed-rounded tip with double lateral vents minimizes the risk of needle binding and wedging in the root canals, trauma to soft/hard tissue and extrusion of the irrigation solution beyond apical foramen

Thin-walled cannula provides improved flow rate for optimal flushing effect

Flexible cannula for increased flexibility and improved access

www.medmix.swiss

medmix US Inc. | 8 Willow Street | Salem, NH 03079, USA

 medmix Deutschland GmbH | Zum Kesselort 53 | 24149 Kiel, Deutschland

Legal Notice: This publication is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. Copyright © 2024 medmix Switzerland AG. All rights reserved.



 **medmix**

DENTAL 