



Painless Steel™

DENTAL INJECTION NEEDLES



PAINLESS STEEL INJECTION NEEDLES

transcendent™ Painless Steel Injection Needles are designed for effective and safe delivery of anesthetic solution to ensure optimal patient comfort.



REFERENCE	DESCRIPTION	SIZE	COLOR	QTY/BOX
162238	30 G Short	0.3 x 25 mm	blue	100
162237	30 G Extra Short	0.3 x 13 mm	blue	100
162240	27 G Long	0.4 x 35 mm	yellow	100
162239	27 G Short	0.4 x 25 mm	yellow	100
162242	25 G Long	0.5 x 35 mm	red	100
162241	25 G Short	0.5 x 25 mm	red	100

PRODUCT FEATURES AND BENEFITS

Siliconized cannula with 3-edge lancet grinding provides smooth entry which minimizes pain and discomfort to the patient

Hub bevel mark facilitates orientation of the bevel

Threaded hub for easy and safe handling

www.medmix.swiss

medmix US Inc. | 8 Willow Street | Salem, NH 03079, USA

 medmix Deutschland GmbH | Zum Kesselort 53 | 24149 Kiel, Deutschland

Legal Notice: This publication is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. Copyright © 2024 medmix Switzerland AG. All rights reserved.



DENTAL 