

MIXPAC™



T-Mixer™ Colibri™

T-Mixer Colibri Plus Breakable

for 18 to 75 mL Systems

- impression of root canal
- impressions of preparations
- core build-ups
- cementation of root post



Ratio	Description	Outer Ø (mm)	Inner Ø (mm)	Gauge	Elements	Ø (mm)	L (mm)*	Vol. (mL)	Order No.	Article
1:1, 2:1	mint, natural elements	1.1	0.8	19	12	4.2	63.8+15	0.62	165 788	MBT 4.2-12-N11B

* length mixing tip + length cannula

T-Mixer Colibri

for 2.5 to 10 mL Systems

- reconstruction with root post
- core build-ups
- cementation of root post and post core
- bulk filling
- base/liner



Ratio	Description	Outer Ø (mm)	Inner Ø (mm)	Gauge	Elements	Ø (mm)	L (mm)*	Vol. (mL)	Order No.	Article
1:1, 2:1	brown, natural elements	0.9	0.6	20	10	2.5	33.85 +15	0.12	154 315	MLT 2.5-10-N09
1:1, 2:1	brown, natural elements	1.1	0.8	19	10	2.5	33.85 +15	0.12	154 316	MLT 2.5-10-N11
1:1, 2:1	brown, natural elements	1.4	1.1	17	10	2.5	33.85 +20	0.14	154 702	MLT 2.5-10-N14

* length mixing tip + length cannula

** Customer specific printing on packaging possible (Standard: Article-No., LOT-No., Made in Switzerland)

Product Advantages

Integrated smoothly rounded stainless steel cannula	<ul style="list-style-type: none"> • Bendable and 360° rotatable - easy access in difficult-to-reach areas. • Constant flow of material, even when the cannula is bent. • No additional intra oral tip is required. • Low risk of soft tissue trauma due to smooth, rounded cannula.
T-Mixer technology	<ul style="list-style-type: none"> • Up to 40% less material waste. • Homogenous mixing results.

medmix Switzerland AG | Ruetistrasse 7 | 9469 Haag, Switzerland
 P +41 81 414 70 00 | www.medmix.swiss



Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC, T-Mixer, and Colibri are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2023



DENTAL