

Protective Sleeves

FOR MIXPAC™ APPLICATION SYSTEMS



MIXPAC Protective Sleeves provide convenient and effective barrier protection of multi-use application systems to protect the dentist, staff and patient against infectious disease during clinical use.

Order No.	Article	Description	Volume	Color	QTY/Box
168058	PRS 801-03	Protective Sleeve, 50 mL, LDPE	50 mL	Transparent	300
168057	PRS 701-03	Protective Sleeve, 5 mL, LDPE	5 mL	Transparent	300
168056	PRS 601-03	Protective Sleeve, 1 mL, LDPE	1 mL	Transparent	300

Product Advantages	Customer Benefits
SAFE	
Disposable single-use barrier	<ul style="list-style-type: none"> • Reduces exposure to chemicals used for cleaning and disinfection. • Reduces the need for cleaning and disinfection.
Tried and tested material	<ul style="list-style-type: none"> • Extensive testing was performed to ensure protection against infectious disease for the dentist, staff and patient. <ul style="list-style-type: none"> - Viral penetration (according to ASTM F 1671) - Synthetic blood penetration (according to ASTM F 1670) - Resistance to puncture (according to ASTM F 1342) - Tear resistance (according to ASTM D 1004) - Tensile strength (according to ASTM D 882) • Biocompatibility: <ul style="list-style-type: none"> - In vitro cytotoxicity (according to EN ISO 10993-5) - Irritation and skin sensitization (according to EN ISO 10993-10) - Systemic toxicity (according to EN ISO 10993-11)
EFFICIENT	
Customized ergonomic design	<ul style="list-style-type: none"> • Snug, conformed fit minimizes visual impairment during clinical use. • Available for 50 mL, 5 mL, 1 mL MIXPAC cartridge systems.
Intuitive and user-friendly	<ul style="list-style-type: none"> • Easily slides over syringe / plunger or cartridge system. • Pre-cut opening easily slips through the mixing tip.
Transparent material	<ul style="list-style-type: none"> • See-through plastic ensures that content inside cartridge is still visible.



medmix USA Inc. | 8 Willow Street | Salem, NH 03079, USA | www.medmix.swiss

Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC Protective Sleeves are manufactured by a qualified medical device manufacturer and distributed by medmix Switzerland AG. MIXPAC is a trademark of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2023



DENTAL