


# Endo Evacuation Needles

## LUER LOCK



### ENDO EVACUATION NEEDLES


transcodent™ Endo Evacuation Needles are designed to reduce trauma to the apex while providing effective suctioning in difficult access areas. The thin wall needle is remarkably flexible, ensuring efficient irrigation even in curved root canals.



REFERENCE	DESCRIPTION	SIZE	COLOR	QTY/BOX
162054	27 G, sterile	0.4 x 32 mm	grey	25
162055	30 G, sterile	0.3 x 32 mm	yellow	25

#### PRODUCT FEATURES AND BENEFITS

- **Apical side vent** for suctioning near the apex
- **4 laser markings** to reduce trauma to the apex
- **Luer lock hub** for use with the Luer Lock Handpieces



medmix US Inc. | 8 Willow Street | Salem, NH 03079, USA | [www.medmix.swiss](http://www.medmix.swiss)

Legal Notice: This publication is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice.  
Copyright © 2023 medmix Switzerland AG. All rights reserved.



DENTAL

