

MIXPAC™

T-Mixer™

EASY TO SAVE MATERIAL



T-Mixer for 18 to 75 mL Systems

- Impression and relining materials
- Bite registration material
- Temporary crown and bridge material
- Restorative materials



Outlet	Ratio	Description	MOQ Box	Inner Ø (mm)	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Article
S	1:1, 2:1	light blue, natural elements	9'000	3.2	12	49.8	0.36	154 610	MBT 3.2-12-S
	4:1,10:1	light blue, orange elements	8'000	3.2	16	57.6	0.46	154 609	MBTX 3.2-16-S
D	1:1, 2:1	light blue, natural elements	9'000	3.2	12	42.5	0.36	128 309	MBT 3.2-12-D
	1:1, 2:1	yellow, natural elements	8'000	4.2	12	50.1	0.59	154 596	MBT 4.2-12-D
	1:1, 2:1	pink, natural elements	6'000	5.4	12	59.4	1.02	154 545	MBT 5.4-12-D
	1:1, 2:1	teal, natural elements	5'000	6.5	12	67.6	1.60	154 568	MBT 6.5-12-D

Product Advantages	Customer Benefits
Innovative mixing tip design allows more compact shape	Reduced material waste volume
Advanced technology leads to more efficient material processing	More homogeneous mixing results enhances application safety
Fully compatible to MIXPAC Systems	No change in working environment

EASY TO SAVE MATERIAL



T-Mixer for 2.5 to 10 mL Systems

- Restorative materials



Outlet	Ratio	Description	Inner Ø (mm)	MOQ Box	Elements	Length (mm)	Waste Vol. (mL)	Order No.	Article
S	1:1, 2:1	brown, natural elements	2.5	9'000	10	34.4	0.13	154 640	MLT 2.5-10-S
	1:1, 2:1	brown/black, natural elements	2.5	9'000	10	34.4	0.13	154 642	MLT 2.5-10-SB
	1:1, 2:1	blue, natural elements	2.5	9'000	10	34.4	0.13	154 622	MLT 2.5-10-SW
	1:1, 2:1	black, natural element	2.5	9'000	10	34.4	0.13	160 646	MLT 2.5-10-S-B
	4:1, 10:1	brown, orange elements	2.5	9'000	10	34.4	0.15	154 625	MLTX 2.5-10-S
	4:1, 10:1	white, orange elements	2.5	9'000	10	34.4	0.15	154 626	MLTX 2.5-10-S-W
	4:1, 10:1	blue, orange elements	2.5	9'000	10	34.4	0.15	154 627	MLTX 2.5-10-SW
D	1:1, 2:1	brown, natural elements	2.5	9'000	10	30.6	0.14	154 636	MLT 2.5-10-D
	1:1, 2:1	purple, natural elements	2.5	9'000	10	30.6	0.14	154 722	MLT 2.5-10-DV
	1:1, 2:1	blue, natural elements	2.5	9'000	10	30.6	0.14	154 637	MLT 2.5-10-DW
	1:1, 2:1	black, natural elements	2.5	9'000	10	30.6	0.14	153 290	MLT 2.5-10-D-B
	4:1, 10:1	brown, orange elements	2.5	9'000	10	30.6	0.18	154 624	MLTX 2.5-10-D

Product Advantages	Customer Benefits
Innovative mixing tip design allows more compact shape	Reduced material waste volume
Advanced technology leads to more efficient material processing	More homogeneous mixing results enhances application safety
Fully compatible to MIXPAC Systems	No change in working environment

medmix Switzerland AG | Ruetistrasse 7 | 9469 Haag, Switzerland
 P +41 81 414 70 00 | www.medmix.swiss



Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC, T-Mixer, and Colibri are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2023



DENTAL