

MIXPAC™

Dispenser

for 5 to 10 mL Systems



5 mL System



Ratio 1:1

Ratio 4:1

10 mL System



Ratio 1:1

Ratio 4:1

Ratio 10:1

MIXPAC Dispensers available for the 5mL and 10mL cartridge systems, help facilitate dispensing and application of high viscosity dental materials by lowering the extrusion force compared to standard plunger and cartridge application systems.

Its ergonomic design provides comfortable, fatigue-free dispensing for precision intraoral application.

5mL System

Item	Ratio	Description	Volume	Color	Order No.	Article
Dispenser	1:1	5 mL Dispenser, including 1:1 Plunger	5 mL	White	112 570	DL 05-01-00
	4:1	5 mL Dispenser, including 4:1 Plunger	5 mL	White	112 569	DL 05-04-00
Plunger	1:1	1:1 Plunger for 5mL Dispenser	5 mL	Black	112 512	PLL 05-01
	4:1	4:1 Plunger for 5mL Dispenser	5 mL	Black	112 513	PLL 05-04

10mL System

Item	Ratio	Description	Volume	Color	Order No.	Article
Dispenser	1:1	10 mL Dispenser, including 1:1 Plunger	10 mL	White	112 568	DL 10-01-00
	4:1	10 mL Dispenser, including 4:1 Plunger	10 mL	White	112 567	DL 10-04-00
	10:1	10 mL Dispenser, including 10:1 Plunger	10 mL	White	112 565	DL 10-10-00
Plunger	1:1	1:1 Plunger for 10mL Dispenser	10 mL	Black	112 514	PLL 10-01
	4:1	4:1 Plunger for 10mL Dispenser	10 mL	Black	112 515	PLL 10-04
	10:1	10:1 Plunger for 10mL Dispenser	10 mL	Black	112 516	PLL 10-10

Product Advantages	Customer Benefits
Reduced extrusion force	• Easy and comfortable application of multi-component, high viscosity dental materials (e.g. core build up materials)
Ergonomic and intuitive design	• Simple, user-friendly assembly • Ideal hand grip form for smaller hands
Robust and durable material	• Strong mechanical performance after sterilization using an autoclave (134°C, 273°F) • High-performance thermoplastic PPS (Polyphenylene Sulfide)
Customization	• Logo printing options available
Harmonized compatibility	• Reliable connection and interaction between all MIXPAC 5mL and 10mL application systems

medmix Switzerland AG | Ruetistrasse 7 | 9469 Haag, Switzerland
P +41 81 414 70 00 | www.medmix.swiss/dental

Legal Notice: This brochure is a general presentation. The information contained in this publication is believed to be accurate and reliable but is not to be construed as implying any warranty or guarantee of performance. We waive any liability and indemnity for effects resulting from its application. Please, contact us for a description of the warranties and guarantees offered with our products. Directions for use and safety will be given separately. All information herein is subject to change without notice. MIXPAC and T-Mixer are trademarks of medmix Switzerland AG. medmix is a trademark of medmix AG. © medmix Switzerland AG 2023

**DENTAL**